

Clinical Bulletin August 2020

Dear Valued Customer,

In our ongoing efforts to further improve our systems we are pleased to announce the release of a new software upgrade for Morpheus8 Applicator.

We are committed to your success and satisfaction and, therefore, it is important for us to ensure that feedback from our customers is implemented to improve the Morpheus8 Handpiece and your user experience.

**This software upgrade includes the following improvement features:**

With this upgrade we are excited to introduce the new Prime Tip (12 pin tip). With a smaller lattice tip of 12 thin gold-coated micropins, and a treatment depth of up to 4mm + an additional heat profile of 1mm, Morpheus8 Prime is the ideal solution for targeting small and hard-to-reach tissue areas, specifically small curved or bony areas such as the periorbital area, forehead, nose, upper lip, inner arms, and other delicate regions. When using the Prime Tip, choose “12 pin” on the screen.



More tips are in R&D pipeline for Morpheus8. For the most updated Quick Reference Guide and Operator Manual please visit our InMode Resource Center at <http://inmoderesources.com/>

We are committed to providing you with the latest innovations in the industry and we can't wait to hear what you think. If you have any questions or concerns, please do not hesitate contact us

Sincerely,

The InMode Clinical Team